

25 years of Renishaw




In 1998, Renishaw celebrates twenty five years of trading, the company having been established by David McMurtry, Chairman and Chief Executive, and John Deer, Deputy Chairman, on 4th April 1973. Today, Renishaw stands at the forefront of automated metrology.

The Group's products provide manufacturers with the ability to machine components accurately, and perform measurement traceable to International Standards. Wherever precision manufacturing is required, the use of Renishaw's products ensures that exacting specifications are met, with cost-effective methods.

Renishaw plc is located at Wotton-under-Edge, Gloucestershire and has established subsidiary companies to provide customer support in its major markets of USA, Japan, Germany, France, Italy, Spain, Switzerland, Hong Kong and Brazil, together with Representative Offices in Singapore, The Peoples' Republic of China, Indonesia and India, and a network of distributors worldwide.

For the financial year ending 30th June 1997, sales totalled £M 81.4, of which £M 72.1 represented business overseas. At that date Renishaw employed 801 in the UK, together with 308 overseas, totalling 1109.

From its leading market position, the Renishaw Group continues to expand its product range into ever increasing business sectors, worldwide. Identifying and targeting new market opportunities has led to the continuous development and introduction of new, highly innovative products which significantly enhance the manufacturing capabilities in a wide range of industries.




Celebrations of Renishaw's Jubilee will be centred on the official opening of the Group's new facilities at New Mills, Wotton-under-Edge in September, which will be followed by a series of "Open Days" for customers. Details will be published nearer the time.

A "Jubilee" Group Profile is available from the company detailing its history as well as current product applications around the globe. Please contact the Spectroscopy Group, or visit our website <http://www.renishaw.com>, if you would like a copy.

continued on back page...

In situ Raman measurements on catalysts

Raman spectroscopy is increasingly being used in the study of catalysts and catalytic reactions. Professor Graeme Hutchings' group at the University of Wales in Cardiff is looking at vanadium phosphorus mixed oxide catalysts (VPOs) used for the selective oxidation of alkanes. Raman spectroscopy is particularly powerful when used in such studies as it can discriminate between the many oxide phases which can exist in the catalysts. In addition, it allows measurements to be done *in situ*, in reactive environments, at high temperature, and because water is relatively weak Raman scatterer, aqueous studies can be undertaken without interference from the water bands.



Raman spectra of various phases of VPO catalyst

continued on page 3...

Ultra-high resolution mapping: a new advance

We are pleased to announce the introduction of a new automated mapping stage for System 1000/2000/3000 Raman microscopes.

The new stage uses another Renishaw product, the RG2 non-contact optical encoder as part of closed-loop feedback systems for the x- and y-axes. This gives unparalleled repeatability and a step size of only 0.1 µm. In other respects the stage performs like the standard 1 µm resolution stage; it is fully integrated within WiRE (Windows-based Raman Environment), Renishaw's proprietary spectroscopy control and analysis software, and can be retrofitted to any Renishaw Raman microscope.

continued on next page...

This edition also contains articles on:

- dual laser solutions
- encoded diffraction grating stage
- long standoff fiber optic probes
- Renishaw Raman at the Laboratory of the Government Chemist

Ultra-high resolution mapping: a new advance *(continued from front cover)*



High resolution motorized mapping stage.
One of the two Renishaw position encoders is circled.


Specifications*

Range of travel:	108 mm x 71 mm, others available
Microscopy:	Leica DMLM or Olympus BHSM
Repeatability:	$\pm 0.5 \mu\text{m}$ ($\pm 5 \mu\text{m}$)
Step size:	$0.1 \mu\text{m}$ ($1 \mu\text{m}$)
Overall dimension of stage:	306 x 225 x 25 mm
Weight of stage:	2.7 kg


*those for the standard (open-loop) stage are given in brackets

The encoded stage allows users to map finer features on their samples and to reliably return to the same position on the sample, even when the stage has been moved appreciable distances from that point.

A typical application that will benefit from the use of the new stage is the mapping of stress round small features on large wafers of integrated semiconductor circuits.



RG2 encoders




WiRE Raman mapping software display reveals the variation in carbon species in a SiC fiber-reinforced composite material

Data courtesy of A. Vassel and M. Parlier (Office National D'Etudes et de Recherches Aeronautiques) and M. Colombari (Centre National de la Recherche Scientifique), France. For further details of this work, please contact the spectroscopy group and request application note AN/065.

Renishaw dual laser solutions

One of the benefits of the high optical efficiency of the Renishaw Raman Microscope is the corresponding reduction in laser power requirements. Small air-cooled lasers can be used, instead of the large water-cooled lasers that were standard a few years ago. Air-cooled lasers are far less expensive (both to buy, and to run), and spectroscopists can now afford two lasers. For example, many Renishaw customers use a 514 nm argon-ion laser for routine work, and a 782 nm laser diode for studying highly fluorescent samples (such as paints).



Easy Change filter mount

Last year Renishaw introduced Easy Change mounts to facilitate the swap of holographic filters that is necessary when the excitation wavelength is changed. These mounts hold the two holographic filters in the correct geometry even when the mount is removed from the instrument.

We now introduce the Dual Laser Mount. This comprises hardware and optics to mount two lasers on the back of a System 1000/2000/3000 instrument. The optics are designed so that the laser beams are colinear when they exit the mount. When you want to change excitation wavelengths you do not have to adjust any hardware on the mount - just turn one laser



Dual Laser Mount with Renishaw 782 nm diode laser (far) and Renishaw 633 nm HeNe laser (near)

off, and the other one on! Together, the laser mounts and the Easy Change filter mounts allow the Renishaw instrument's laser wavelengths to be swapped in a couple of minutes.

The Dual Laser Mounts are available for a range of lasers; please contact your local representative for more details.


***In situ* Raman measurements on catalysts**

(continued from front cover)

Oxide catalysts are used in a broad range of important industrial processes and are generally prepared in the form of a precursor which is converted into the active catalyst under well defined reaction conditions. VPO catalysts are used to effect the oxidation of n-butane to maleic anhydride. $(VO)_2P_2O_7$ is considered to be the main active component, and is prepared from the precursor $VOHPO_4 \cdot 0.5H_2O$. Hutching's group is studying the conversion of the precursor to the active form in real time and under the true reaction conditions. Characterisation of the catalyst under a flow of butane in air has been achieved between room temperature and 1000 °C, showing the interconversion between the many phases of the VPO catalyst. The gaseous out flow is then passed into a gas chromatograph to analyse the product mixture. Hence a complete description of the catalytic process including the evolution of the active catalyst is obtained. This system has widespread use for many catalytic systems, and further work will involve the addition of promoters and modifiers to the catalyst and the gas stream to study their effects on the formation of the active catalytic species. Further studies are envisaged using the fiber optic probe system to monitor reactions *in situ* at elevated pressure. This should further enhance the group's ability to study catalyst systems under true reaction conditions.

For further details of this work please contact:
Dr Richard Wells, Chemistry Department
University of Wales at Cardiff.
email: WellsRP@cardiff.ac.uk

Renishaw encoded diffraction grating stage



Renishaw encoded diffraction grating stage


The RG2 optical encoder is also used in the Renishaw encoded rotary stage. This stage is used as a mount for the diffraction grating in many Renishaw spectrometers. It gives a far higher repeatability than previously possible, typically 2 arcsec. This corresponds to a spectral repeatability of $\pm 0.25 \text{ cm}^{-1}$ when a 1800 groove mm^{-1} grating is used on a System 1000 Raman microscope.

The encoded diffraction grating stage is fitted as standard to all new spectrometers, and can be retrofitted to existing Renishaw System 1000/2000/3000 instruments at the customer's site.

Long standoff fiber optic probes

Microscope objectives are ideal lenses for Raman spectroscopy in most respects; they give high magnification images, and their high numerical apertures result in high collection efficiencies. However they normally have working distances (distance from the front of the lens to the sample) of only a few mm. This makes them unsuitable for some applications, such as measurements inside reaction chambers and cells.

The Renishaw spectroscopy group has been developing long working distance optics for the Renishaw fiber optic probes. The spectrum below illustrates the high performance achievable.



Raman spectrum of propan-2-ol, exposure time 2 s

Conditions

Spectrometer	Renishaw System 100
Excitation wavelength	514.5 nm (Ar ⁺)
Power on sample	8.5 mW
Working distance	0.55 m
Optical fiber length	10 m
Collection lens aperture	25 mm

Please contact us if such long standoff optics would help you in your applications.

Renishaw Raman at Laboratory of the Government Chemist



Renishaw has recently supplied one of its Raman systems to the Laboratory of the Government Chemist (Teddington, UK) to establish a national (and international) facility for Raman micro-analysis. The system is available to customers requiring an analytical service on either a one off sample/problem-solving basis or on a longer term structured research basis.

For further information please contact either:
Dr. Henryk Herman, LGC, +44 (0)181 943 7316
Dr. Ken Williams, Renishaw plc, +44 (0)1453 844302

25 Years of Renishaw (continued from front cover)

- 1991 Renishaw develops the first prototype of a new generation of compact Raman microscope instrumentation
- 1992 Raman imaging microscope launched
Renishaw sells 1st Raman imaging microscope
- 1993 Class 1 laser-safe version of Raman imaging microscope
- 1994 782 nm diode laser - the first Raman-ready diode laser
- 1995 Raman fibre optic probe for remote analysis
electrochemical cells for *in situ* Raman measurements
- 1996 RenCam, Renishaw's CCD camera
Renishaw encoded grating stage
- 1997 UV Raman microscope, DLC analyzer, encoded mapping stage
- 1998 System 100 - a compact Raman process analyzer,
Raman libraries, Renishaw HeNe laser

Key milestones of the Renishaw spectroscopy group

Meet the people

We plan to introduce some of the spectroscopy group to you in each issue, so that you can put faces to the voices you hear on the telephone. To start, here are some of the UK marketing group personnel:



L to R: Anne Bell, Anette Zimmermann, Ian Hayward, Chris Dyer, Jon Woolls, Ian Wilcock, Ken Williams.

This group supports world-wide sales activities and our subsidiaries (see below) and distributors.

Renishaw worldwide companies with spectroscopy support**Renishaw GmbH**

Karl-Benz Strasse 12,
D 72124 Pliezhausen, Germany
Tel +49 712 797960
Fax +49 712 788237

Renishaw S.p.A.

Via dei Prati 5, 10044 Pianezza,
Torino, Italy
Tel +39 11 966 10 52
Fax +39 11 966 40 83

Renishaw S.A.

15 rue Albert Einstein,
Champs sur Marne,
77437 Marne la Vallee, Cedex 2,
France
Tel +33 1 64 61 84 84
Fax +33 1 64 61 65 26

Renishaw Inc.

623 Cooper Court,
Schaumburg, Illinois 60173, USA
Tel +1 847 843 3666
Fax +1 847 843 1744

USA - East Coast

Renishaw Inc
1493 Beacon Street
Brookline MA 02146
Tel: 1 617 975 0245
Fax: 1 617 975 0269

USA - West Coast

Renishaw Inc
2992 Neal Avenue
San José CA 95128
Tel: 1 408 260 9430
Fax: 1 408 260 9431

Renishaw (Hong Kong) Ltd

Unit 4A, 3/F, New Bright Building,
11 Sheung Yuet Road,
Kowloon Bay, Hong Kong
Tel +852 2753 0638
Fax +852 2756 8786

The Peoples' Republic of China,

Room 2302 Silver Tower,
No 2 Dong San Huan North Road,
Beijing 100027
Tel +86 10 6410 7993
Fax +86 10 6410 7792

If you or any of your colleagues would like to be added to our mailing list or if you have any queries or comments about the contents of this newsletter, please contact us at any of the addresses above or simply email us at raman@renishaw.co.uk.

Your local Renishaw Spectroscopy agent or distributor is:

Renishaw plc, Transducer Systems Division,
Old Town, Wotton-under-Edge,
Gloucestershire GL12 7DH, United Kingdom
Tel +44 (0) 1453 844302 Fax +44 (0) 1453 844236
email: spectroscopy.products@renishaw.com
<http://www.renishaw.com>