

## 使用 EDGE 加压流体萃取仪快速提取不同样品中的 PFAS

PFAS (全氟和多氟烷基物质) 是一种广泛应用于工业生产中的人造化合物，其耐高温，防污防滑的特性使得它在不粘炊具，食品快餐盒，消防泡沫，防污织物和家具中被广泛应用。因极难在环境中降解，具有长距离迁移及生物累积性的特点，使得 PFAS 一旦进入人体，就会一直留在人体之中。PFAS 类物质的累积会引起诸多人类健康问题，如癌症、内分泌紊乱和不孕症等。因此，对他们的监测至关重要。

最近 CEM 北美实验室用 EDGE 加压流体萃取系统分别对随机购买的黄瓜、纸杯蛋糕及微波食品进行 PFAS 的提取，同时以土壤为样本进行低中高浓度的 PFAS 加标实验。

为确保实验不受 PFAS 干扰，已对 EDGE 系统中可能存在 PFAS 的材料进行了全面替换，并在测试过程中，避免使用可能含有 PFAS 的耗材及检测设备。

**Table 1.** PFAS Compounds Assessed, Their Respective Detection Limit, and Soak Test Results

Compound	Limit of Detection (ng/L)	EDGE System	PEEK Tubing	Polypropylene Tubing	C9 Q-Disc	G1 Q-Disc	Q-Disc Separator	Q-Matrix Hydra	Q-Cup
9-chlorohexadecafluoro-3-oxanone-1-sulfonic acid (9Cl-PF30NS)	25	ND	ND	ND	ND	ND	ND	ND	ND
11-chloroeicosafluoro-3-oxaundecane-1-sulfonic acid (11Cl-PF30UdS)	25	ND	ND	ND	ND	ND	ND	ND	ND
1H, 1H, 2H, 2H-perfluorodecane sulfonic acid (8:2 FTS)	50	ND	ND	ND	ND	ND	ND	ND	ND
1H, 1H, 2H, 2H-perfluorooctane sulfonic acid (6:2 FTS)	50	ND	ND	ND	ND	ND	ND	ND	ND
1H,1H,2H,2H-perfluorododecane sulfonic acid (10:2 FTS)	50	ND	ND	ND	ND	ND	ND	ND	ND
1H,1H,2H,2H-perfluorohexane sulfonic acid (4:2 FTS)	50	ND	ND	ND	ND	ND	ND	ND	ND
Hexafluoropropylene oxide dimer acid (GenX)	50	ND	ND	ND	ND	ND	ND	ND	ND
4,8-dioxa-3H-perfluorononanoic acid (ADONA)	50	ND	ND	ND	ND	ND	ND	ND	ND
N-ethylperfluoro-1-octanesulfonamide (EtFOSA)	50	ND	ND	ND	ND	ND	ND	ND	ND
N-ethylperfluoro-1-octanesulfonamidoacetic acid (EtFOSAA)	50	ND	ND	ND	ND	ND	ND	ND	ND
2-N-ethylperfluoro-1-octanesulfonamido-ethanol (EtFOSE)	50	ND	ND	ND	ND	ND	ND	ND	ND
N-methylperfluoro-1-octanesulfonamide (MeFOSA)	100	ND	ND	ND	ND	ND	ND	ND	ND
N-methylperfluoro-1-octanesulfonamidoacetic acid (MeFOSAA)	50	ND	ND	ND	ND	ND	ND	ND	ND
2-N-methylperfluoro-1-octanesulfonamido-ethanol (MeFOSE)	50	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-1-butanedisulfonic acid (PFBS)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-1-decanedisulfonic acid (PFDS)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-1-heptanedisulfonic acid (PFHpS)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-1-nonanedisulfonic acid (PFNS)	53	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-1-octanesulfonamide (PFOSA)	50	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-1-pentanesulfonic acid (PFPeS)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluorohexanesulfonic acid (PFHxS)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-butanoic acid (PFBA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-decanoic acid (PFDA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-dodecanoic acid (PFDoA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-heptanoic acid (PFHpA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-hexanoic acid (PFHxA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-nonanoic acid (PFNA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-octanoic acid (PFOA)	25	ND	ND	ND	ND	ND	ND	ND	ND
Perfluoro-n-pentanoic acid (PFPeA)	25	ND	ND	ND	ND	ND	ND	ND	ND



### PFAS Detected

2-N-methylperfluoro-1-octanesulfonamido-ethanol  
(MeFOSE)

perfluoro-n-butanoic acid (PFBA)

perfluoro-n-hexanoic acid (PFHxA)

perfluoro-n-octanoic acid (PFOA)

perfluoro-n-pentanoic acid (PFPeA)



### PFAS Detected

N-ethylperfluoro-1-octanesulfonamidoacetic acid  
(EtFOSAA)

N-methylperfluoro-1-octanesulfonamidoacetic acid  
(MeFOSAA)

perfluoro-1-pentanesulfonic acid (PFPeS)

perfluorohexanesulfonic acid (PFHxS)

perfluoro-n-butanoic acid (PFBA)

perfluoro-n-heptanoic acid (PFHpA)

perfluoro-n-hexanoic acid (PFHxA)

perfluoro-n-nonanoic acid (PFNA)

perfluoro-n-octanoic acid (PFOA)

perfluoro-n-pentanoic acid (PFPeA)



PFAS Detected	
N-ethylperfluoro-1-octanesulfonamidoacetic acid (EtFOSAA)	perfluoro-n-hexanoic acid (PFHxA)
N-methylperfluoro-1-octanesulfonamidoacetic acid (MeFOSAA)	perfluoro-n-nonanoic acid (PFNA)
perfluoro-1-butanefulfonic acid (PFBS)	perfluoro-n-octanoic acid (PFOA)
perfluoro-1-decanesulfonic acid (PFDS)	perfluoro-n-pentanoic acid (PFPeA)
perfluoro-1-heptanesulfonic acid (PFHpS)	perfluoro-n-butanoic acid (PFBA)
perfluoro-1-pentanesulfonic acid (PFPeS)	perfluoro-n-decanoic acid (PFDA)
perfluorohexanesulfonic acid (PFHxS)	perfluoro-n-dodecanoic acid (PFDoA)
perfluoro-n-butanoic acid (PFBA)	perfluoro-n-heptanoic acid (PFHpA)
perfluoro-n-decanoic acid (PFDA)	perfluoro-n-tridecanoic acid (PFTrDA)
perfluoro-n-dodecanoic acid (PFDoA)	perfluoro-n-u0ecanoic acid (PFUdA)
perfluoro-n-heptanoic acid (PFHpA)	perfluorooctanesulfonic acid (PFOS)



Compound	Low Spike	RSD (n=3)	High Spike	RSD (n=3)
1H, 1H, 2H, 2H-perfluorodecane sulfonic acid (8:2 FTS)	83.00%	7.81%	87.00%	7.00%
1H, 1H, 2H, 2H-perfluorooctane sulfonic acid (6:2 FTS)	80.00%	3.00%	87.67%	7.09%
1H,1H,2H,2H-perfluorohexane sulfonic acid (4:2 FTS)	86.67%	10.02%	87.67%	5.13%
N-ethylperfluoro-1-octanesulfonamidoacetic acid (EtFOSAA)	85.67%	6.66%	86.00%	5.57%
N-methylperfluoro-1-octanesulfonamidoacetic acid (MeFOSAA)	78.67%	1.15%	81.33%	4.16%
Perfluoro-1-butanedisulfonic acid (PFBS)	78.00%	3.00%	86.00%	1.73%
Perfluoro-1-decanedisulfonic acid (PFDS)	77.00%	6.56%	84.33%	0.58%
Perfluoro-1-heptanedisulfonic acid (PFHpS)	81.00%	2.65%	87.00%	2.65%
Perfluoro-1-nonanedisulfonic acid (PFNS)	76.00%	6.56%	84.67%	2.89%
Perfluoro-1-octanesulfonamide (PFOSA)	81.00%	8.00%	86.67%	5.03%
Perfluoro-1-pentadisulfonic acid (PFPeS)	78.00%	4.36%	86.33%	0.58%
Perfluorohexanedisulfonic acid (PFHxS)	98.33%	2.89%	89.67%	3.79%
Perfluoro-n-butanoic acid (PFBA)	85.65%	4.01%	88.17%	2.20%
Perfluoro-n-decanoic acid (PFDA)	101.33%	7.51%	93.33%	2.89%
Perfluoro-n-dodecanoic acid (PFDoA)	79.67%	7.37%	79.33%	2.08%
Perfluoro-n-heptanoic acid (PFHpA)	90.67%	9.02%	81.67%	3.21%
Perfluoro-n-hexanoic acid (PFHxA)	100.67%	9.02%	91.08%	6.64%
Perfluoro-n-nonanoic acid (PFNA)	96.33%	3.21%	92.67%	2.52%
Perfluoro-n-octanoic acid (PFOA)	82.77%	3.87%	85.53%	2.25%
Perfluoro-n-pentanoic acid (PFPeA)	79.50%	0.71%	85.67%	4.04%
Perfluoro-n-tetradecanoic acid (PFTeDA)	86.67%	4.73%	88.67%	1.15%
Perfluoro-n-tridecanoic acid (PFTrDA)	63.00%	3.46%	68.33%	3.51%
Perfluoro-n-undecanoic acid (PFUdA)	76.33%	1.53%	79.67%	3.51%
Perfluorooctanesulfonic acid (PFOS)	84.00%	4.00%	79.60%	0.53%

EDGE 对低、中、高浓度加标的土壤样品的提取，单样品提取时间小于 10 分钟，并可获得理想回收率及 RSD 值，同时，在黄瓜、纸杯蛋糕、微波食品中均检出 PFAS 物质。