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## Identification Tetramine' Formulas by MassWorks Combination of Unit-resolved Mass Spectrometry

Li Weijian, Gong Junjie, Ouyang Weimin  
(Lumiere Tech Ltd ,Beijing 100080)

**Abstract:** This application example demonstrates that the comprehensive mass spectral calibration involving both mass and peak shape enables formula determination on the GC/MSD system. The MassWorks™ calibration can be conveniently established using the on-board PFTBA tune standard through infusion measurement within the same GC/MS data acquisition. While the mass accuracy achievable on the system may not be the industry standard 5ppm for all ions at different m/z values, the combination of mass accuracy and Spectral Accuracy can indeed achieve formula determination on a real chromatographic time scale. Not only can this formula ID capability be applied to unknown identifications where there is no available library spectrum to compare with through conventional library search, it can also be utilized as an independent confirmation of library search results through formula ID of either the molecular ion or its fragment ions.

**Key words:** MassWorks™, accurate mass, mass spectral calibration ,Formula determination, unit mass resolution.