## APPLICATION NOTE

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DHA Applications	
<div>  New            Edit            Reimport            Delete            Export            Close         </div>	
Application Name	Description
▶ ASTM D 7900	<p>This test method specifies a method to determine the boiling range distribution of hydrocarbons in stabilized crude oil up to and including n-nonane. A stabilized crude oil is defined as having a Reid Vapor Pressure equivalent to or less than 82.7 kPa. The results of this test method can be combined with those from Test Method D7169 and IP 545 to give a full boiling point distribution of a crude oil. RON and MON values are not calculated and set to 0.</p> <p>Oven 35°C (30 min) =&gt; 2 °C/min =&gt; 200 °C (20 min)            PTV 35°C (0.1 min) =&gt; 700 °C/min =&gt; 200 (30 min) =&gt; (700) =&gt; 250 °C            Total Flow 93.9 ml/min Split Ratio 1:100 Split Flow 90 ml/min Helium            Inj Vol 1.0 µl (5 µl syringe)            Col1: 5 m 0.25mm 1µm polydimethylsiloxane            Col2: 50m 0.21mm 0.5µm polydimethylsiloxane            Col1: Pressure 37.7 psi Flow 1.2 [ml/min] (40 min) =&gt; 99 ml/min =&gt; -20 ml/min            Col2: Pressure 35.8 psi Const Flow 1.5 [ml/min]</p>
ASTM D 6730	<p>This test method covers the determination of individual hydrocarbon components of spark-ignition engine fuels and their mixtures containing oxygenate blends (MTBE, ETBE, ethanol, and so forth) with boiling ranges up to 225°C. Other light liquid hydrocarbon mixtures typically encountered in petroleum refining operations, such as blending stocks (naphthas, reformates, alkylates, and so forth) may also be analyzed. Individual component concentrations must be in the range from 0.01 to approximately 30 mass %.</p>
ASTM D 6729	<p>This test method covers the determination of individual hydrocarbon components of spark-ignition engine fuels and their mixtures containing oxygenate blends (MTBE, ETBE, ethanol, and so forth) with boiling ranges up to 225°C.</p>

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