


# Lenses or No Lenses? A Study of Ion Transfer Efficiency at Interfaces in a Lens-free Triple Quadrupole MS



ASMS 2011, TP04

**Felician Muntean and Roy Moeller**  
Bruker Chemical and Applied Markets  
3500 West Warren Avenue  
Fremont, CA 94538 USA

## Introduction

The efficiency of ion transfer between the quad mass analyzers and the  $Q_2$  collision cell in a triple quad is critical for instrument robustness and sensitivity; qualities that make triple quads "work-horses" in quantitative MS today.


### Challenging requirements

- $Q_2$  contains collision gas but there should be no collisions in  $Q_1$  or  $Q_3$ , which are in close proximity, in the same vacuum region.
- The ion beam is wide/divergent at the output of each quad but the next quad requires a narrow, collimated beam.

### Common implementations today

- Traditional - Lens approach: collision cell around  $Q_2$  with narrow apertures in and out to limit gas flow. Electrostatic lenses used to refocus ions in and out. Common but requires complex tuning, prone to contamination, nodding.
- Newer - Lens-free approach (Fig. 5): collision area is defined with seals around and outside of  $Q_2$  ion guide path (US Pat 6576897). RF-only ion guides and not lenses are used to guide ions in and out of  $Q_2$ . Can provide high transmission, simple/robust instrument design.

This computational work studies the ion transfer efficiency and fundamental limitations of the two approaches above.



## Methods

### Quad interface models (Fig. 1)

- Lens-free ( $Q_1Q_2$  and  $Q_2Q_3$ ), Fig. 1a. Gas cell boundaries are marked with dashed lines. Ion transfer is achieved through close proximity RF-only ion guides:  $Q_1$  post-filters,  $Q_2$  RF-only input and output sections and  $Q_3$  pre-filters.
- Three-lenses ( $Q_1LLLQ_2$  and  $Q_2LLLQ_3$ ), Fig. 1b. Gas is contained in  $Q_2$  with one narrow aperture and 2 other ion lenses at each end for refocusing the ion beam.  $Q_1$  post-filters and  $Q_3$  pre-filters.

Ion trajectory simulations are done using SIMION 8 (SIS Inc., Ringoes, NJ 08551; www.simion.com) and LUA programming.

#### Simulation Details

Geometry (Fig. 1): Round rod quads  $r_0$  4 mm and length 40 mm; pre-post filter length 12.5 mm; square-rod  $Q_2$   $r_0$  3 mm, sides 4.5 mm, length 40 mm. 1.5 mm per grid point. Axial distance between all elements 1 mm. For  $Q_1LLLQ_2$ , lens 1 ID 6 mm, length 4 mm, lens 2 ID 6 and 4 mm, length 15 and 1 mm, lens 3 ID 2 mm, length 2 mm. For  $Q_2LLLQ_3$ , lens 1 ID 2 mm, length 2 mm, lens 2 ID 5 and 3 mm, length 9 and 1 mm, lens 3 ID 3 mm, length 3 mm.


Voltagess: RF frequency 1 MHz, quad RF voltages 306 V at  $m/z$  264, DC voltages adjusted for unit mass separation; q-values pre/post filter 0.5, q-value  $Q_2$  0.5. Tuning voltages without collisions are in Table 1, with collisions  $Q_2$  axis voltage is added to all following elements.

Ion beam: Circle distribution 1 mm diameter, ion energy 5 eV. Trajectories: 96 ions, 8 RF phases per simulation point. Ion  $m/z$  264 or as described in the text.

Collisions: Elastic, hard sphere model, collision gas Ar (40 Da), pressure 0.266 Pa. Collision cross section  $220 \text{ \AA}^2$ . Collisions are assumed to be present outside the collision cell a distance equal to the diameter of collision cell hole.

Configuration	Q1	Q1 post-filter	L1	L2	L3	Q2	L1	L2	L3	Q3	Q3 pre-filter
Q1-Q2	-5	-5				-5					
Q2-Q3						-5				-5	-5
Q1-LLL-Q2	-5	-8	-200	-10	-200	-5					
Q2-LLL-Q3						-5	-100	-10	-75	-20	-5

Table 1. Tuned voltages for quad axes and lenses at  $m/z$  264.



## Results




Fig. 2 Details of ion trajectories in and out of  $Q_2$  for the lens-free (a) and three-lens (b) quad interfaces with and without gas in  $Q_2$ .

### Ion trajectory analysis. Lens-free quad interface (Fig. 2a)

- Smooth transfer between RF-only quad sections.
- No losses due to collisions when  $Q_2$  gas is on. Collisions occur within RF-only guides, they help collimate the ion beam and improve transmission efficiency.
- Most ion losses occur in  $Q_3$ , due to limited mass filter acceptance.
- No significant electrostatic focusing can be achieved with the pre-post filter ion guides.

### Three-lens quad interface (Fig. 2b)

- Efficient ion focusing and transfer when  $Q_2$  gas is off and no collisions occur.
- Significant losses when collision gas is on. Most ion losses occur within the electrostatic lenses, especially in the aperture region.
- Lenses cannot efficiently focus ions in a collision environment.

### Ion transmission analysis (Fig. 3)

Voltages in both systems are "tuned" vs. mass for maximum efficiency in absence of collisions. When collision are on, voltages following  $Q_2$  are adjusted to account for ion energy loss.

- The ion transfer out of  $Q_2$  is less efficient than into  $Q_2$  due to the smaller acceptance of analyzing quad  $Q_3$  compared to the RF-only quad  $Q_2$  (see Fig. 3,  $Q_2 \dots Q_3$  vs.  $Q_1 \dots Q_2$ ).
- The lens-free transmission is rather unaffected by collision pressure, energy or mass, except at low masses  $< m/z$  200.
- The 3-lens transmission is significantly affected by collision pressure, energy, and mass.
- The 3-lens transmission improves with collision energy due to reduced scattering losses.




Fig. 3 Transmission efficiency in and out of  $Q_2$  for the lens-free and the 3-lens interfaces as a function of  $Q_2$  pressure (collision energy 25eV and  $m/z$  264), collision energy (pressure 0.266Pa and  $m/z$  264) and mass (collision energy 25eV and pressure 0.266Pa).

### Q3 low res.




Fig. 4 Transmission efficiency  $Q_2 \dots Q_3$  when  $Q_3$  is set to lower resolution.

- As collision pressure increases, the lens-free interface is significantly more efficient than the 3-lens interface (Fig. 3a).
- When  $Q_3$  acceptance is increased by reducing  $Q_3$  resolution, the lens-free interface provides mostly 100% efficiency but the efficiency of the 3-lens interface is still limited due to collision losses (Fig. 4).




Fig. 5 An example - the lens-free  $Q_2$  and  $Q_2Q_3$  coupling in the Bruker 320 triple-quad.

## Conclusions

- Gas collisions may increase the efficiency of a lens-free, RF-only interface between a quad and collision cell but they significantly reduce the efficiency of lens-based interfaces.
- Collision pressure, energy, and mass have little influence to the ion transfer efficiency of a lens-free interface while they significantly change the efficiency of lens-based interfaces
- Overall, a lens-free, RF-only interface is not just simpler and more robust than a lens-based interface but it also provides significantly higher efficiency.

Triple Quad MS