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

Figure 1. Evolution of the surface heterogeneity of lactose powders with micronization and conditioning

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- ref 1. E. Brendlé, E. Papirer, Surface Properties Characterization by Inverse Gas Chromatography (IGC) Applications, Powders and Fibers, 2006, 47-122  
ref 2. Adscientis, Technical sheet on IGC-ID 2017  
ref 3. Adscientis, Technical sheet on IGC-FC 2017  
ref 4. M.D. Ticehurst, P. York, R.C. Rowe, S.K. Dwivedi, Characterisation of the surface properties of a lactose monohydrate with inverse gas chromatography, used to detect batch variation, International Journal of Pharmaceutics 141, 1996, 93-99  
ref 5. H.E. Newell, G. Buckton, D.A. Butler, F. Thielmann, D.R. Williams, The use of inverse phase gas chromatography to measure the surface energy of crystalline, amorphous and recently milled lactose, Pharmaceutical Research 18, 2001, 662-666  
ref 6. H.E. Newell, G. Buckton, D.A. Butler, F. Thielmann, D.R. Williams, The use of inverse phase gas chromatography to study the change of surface energy of amorphous lactose as a function of relative humidity and the processes of collapse and crystallization, International Journal of Pharmaceutics 217 2001, 45-56  
ref 7. H. Balard, Estimation of the surface energetic heterogeneity of a solid by inverse gas chromatography, Langmuir, 13 1997, 1260-1270.