

Applications of Mixture Analysis to FT-IR Spectra of Multi-Component Polymeric Systems



Dana Garcia,¹ and Marie Scandone,²
(1) Arkema Inc. 900 First Avenue, King of Prussia, PA 19406-0936, (2) Bio-Rad Laboratories, Inc., Informatics Division, Two Penn Center Plaza, Suite 800, 1500 John F. Kennedy Blvd., Philadelphia, PA 19102-1737

Introduction

IR spectral identification of multi-component samples is a challenging process:

- Highly overlapping peaks
- Lack of specific component knowledge
- Assignment of functional groups to a specific component
- Time intensive one-component-at-a-time spectrum analysis

KnowItAll® mixture analysis effectively addresses IR spectral identification of multi-component samples.

Mixture analysis is a mathematical process to extract:

- Combinations of IR spectra in a library database
- Simultaneously minimizing the residual relative to the unknown IR spectrum


Mixture analysis results:

- Provide a series of potential unknown component IR spectra for evaluation
- Improve data mining efficiency
- Simplify search results evaluation

Mixture Analysis


Algorithm for Spectral Fitting:

Finding combinations of database IR spectra that minimize residual.




Residual intensity comparable or smaller than obtained through standard spectral subtraction.

Example 1. Two Component Polymeric System




Component DA#9858 could not be identified through simple searching.

Example 2. Polymer Plasticizer Blend




Both components identified in a single step.

Example 3a. Three Component Polymeric Filler System



Example 3b. Three Component Polymeric Filler System Filler Identification



A combination of mixture analysis and simple search identified all three components of the formulation:

- Tecnoflon® & Paracril® 3300 & Dicalite® SA-23

Conclusions

IR spectral identification of multi-component samples lends itself to a combination approach using mixture analysis and simple spectral search.

High degree of module interconnectivity.

Integrated with all other functionalities of KnowItAll®.

- Individual module provide partial answers
- In totality a systematic and efficient analysis

Applicable to a wide variety of spectral identification problems.

Tecnoflon® is a registered trademark of Solvay Solexis, Paracril® is a registered trademark of Paratec Elastomers LLC
Pebax® is a register trademark of Arkema France

Disclaimer: The statements, technical information and recommendations contained herein are believed to be accurate as of the date hereof. Since the conditions and methods of use of the information referred to herein are beyond our control, Arkema expressly disclaims any and all liability as to any results obtained or arising from any reliance on such information; NO WARRANTY OF FITNESS FOR ANY PARTICULAR PURPOSE, WARRANTY OF MERCHANTABILITY, OR ANY OTHER WARRANTY, EXPRESS OR IMPLIED, IS MADE CONCERNING THE INFORMATION PROVIDED HEREIN. The user should thoroughly test any application before commercialization. Nothing contained herein constitutes a license to practice under any patent and it should not be construed as an inducement to infringe any patent, and the user is advised to take appropriate steps to be sure that any proposed action will not result in patent infringement.
© 2010 Arkema Inc.